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#### Non-Dengue Infection for IgM/IgG test results

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- Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12
- Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Microbiol Immunol. 1994; 38(12): 951-5.

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# GB's – RF

(LATEX AGGLUTINATION TEST)

KIT NAME	KIT SIZE	CAT. NO
GB's RF	25 Tests	SRF0000025T
GB's RF	50 Tests	SRF0000050T
GB's RF	100 Tests	SRF0000100T

## INTRODUCTION

Measurement of rheumatoid factor is used for differentiating rheumatoid arthritis from other chronic inflammatory arthritis and is important in the progress and therapeutic management of the disease. Rheumatoid factor has been associated with some bacterial and viral infections (ex. Hepatitis, infectious, Mononucleosis) some chronic infections (ex. Tuberculosis, parasitic disease, sub acute Bacterial Endocardities) and cancer.

## METHOD PRINCIPLE

The latex reagent coated with the Human gammaglobulin (gG). The test specimen (serum) is mixed with RF latex reagent and allowed to react. If RF is present within detectable levels then a visible agglutination is observed. If RF is absent below detectable levels then no agglutination is observed.

## REAGENTS

Reagent Name	SRF0000025T	SRF0000050T	SRF0000100T
R1 RF latex	1 vial	1 vial	1 vial
R2 Positive control	1 vial	1 vial	1 vial
R3 Negative control	1 vial	1 vial	1 vial

## WORKING REAGENT PREPARATION AND STABILITY

1. Store the reagent at 2-8°C. DO NOT FREEZE.
2. The shelf life of the reagent is as per the expiry date mentioned on the reagent vial labels.

## SPECIMEN

Only serum should be used to testing. Should a delay in testing occur, store the samples at 2-8°C. Samples can be stored for upto a week. Do not use hemolysed serum.

## MATERIAL PROVIDED WITH THE KIT

**Reagents:** RF latex reagent, positive control, Negative control  
**Accessories:** Glass Slide, Plastic Droppers, Mixing sticks.

## ADDITIONAL MATERIAL REQUIRED

Stop watch, high intensity to room temperature before testing.

## NOTES:

1. In vitro diagnostic reagent for laboratory and professional use only. Not for medicinal use.
2. All the reagents derived from human source have been tested for HBsAg and Anti HIV antibodies and are found to be non-reactive. However handle the material as if infectious.
3. Reagent contains 0.1% sodium Azide as a preservative. Avoid contact with skin and mucosa. On disposal flush with large quantities of water.
4. The reagent can be damaged due to microbial contamination or on exposure to extreme temperatures. It is recommended that the performance of the reagent be verified with the positive and negative controls supplied with the kit.
5. Shake the RF latex reagent well before use to disperse the latex particles uniformly and improve test readability.

6. Only a clean and dry glass slide must be used. Clean the slide with distilled water and wipe dry.
7. Accessories provided with the kit only must be used for optimum results.

## PROCEDURE

Bring reagent and samples to room temperature before use.

### Qualitative Method

1. Pipette one drop of serum onto the glass slide using the disposable plastic droppers provided with kit.
2. Add one drop of RF latex reagent to the drop of test specimen on the slide. Do not let the dropper tip touch the liquid on the slide.
3. Using a mixing stick, mix the serum and RF latex reagent uniformly over the entire circle.
4. Immediately start a stop watch. Rock the slide gently, back and forth observing for agglutination macroscopically at two minutes.

### Semi Quantitative Method:

1. Using normal saline prepare serial dilutions of the Serum sample positive in the qualitative method 1:2, 1:4, 1:8, 1:16, 1:32, 1:64 and so on.
2. Pipette each dilution of the serum sample onto separate reaction circles.
3. Add one drop of RF latex reagent to each drop of the diluted serum sample on the slide. Do not let the dropper tip touch the liquid on the slide.
4. Using a mixing stick, mix the sample and the latex reagent uniformly over the entire circle.
5. Immediately start a stop watch. Rock the slide gently, back and forth, observing for agglutination macroscopically at two minutes.

## INTERPRETATION OF TEST RESULTS

### QUALITATIVE METHOD

Agglutination is a positive test result and indicates the presence of detectable levels of RF in the test specimen.  
 No agglutination is a negative test result and indicates the absence of detectable levels of RF in the test specimen.

### SEMI QUANTITATIVE METHOD

Agglutination in the highest serum dilution corresponds to the amount of RF in IU/ml present in the test specimen.

### CALCULATION

$$RF \text{ (IU/ml)} = S \times D$$


Where S = Sensitivity of the reagent i.e., 12 IU/ml  
 D = Highest dilution of serum showing agglutination.

**REMARKS:**

1. Marked lipemic, hemolysed and contaminated serum samples could produce non-specific results.
2. Use of plasma rather than serum can lead to false positive results.
3. Do not read results beyond two minutes.
4. Rheumatoid factors are not exclusive found in rheumatoid arthritis but sometimes in syphilis, systemic lupus erythematosus, hepatitis, hepergammaglobulinemia also.
5. It is recommended that results of the test should be correlated with clinical findings to arrive at the final diagnosis.
6. RF reagent is sensitive to the presence of IgM RF with heterogeneous specificity.

**LITERATURE**

1. Heller G., Jacobson S.A., Koloday MH, Kammerer M.H., J. Immunol, 72:66(1954).
2. Singer J.M, bull Rheum, Dis, 24: 762 (1974).

**Genuine Biosystem Private Limited**

Plot No.97 & 98, kattabomman street,  
Parvathy Nagar Extension,  
Old Perungalathur, Chennai - 600063, India.  
Ph: +91-44-48681845

# GB's - ASO (LATEX AGGLUTINATION TEST)

KIT NAME	KIT SIZE	CAT. NO
GB's ASO	25 Tests	SASO000025T
GB's ASO	50 Tests	SASO000050T
GB's ASO	100 Tests	SASO000100T

## INTRODUCTION

The first group of streptococci are called hemolytic streptococci which can be further subdivided into group (a), group (b), group (c) and group (d). It includes most of the species associated with primary streptococcal infections in humans. The group (a) hemolytic streptococci produce various exotoxins such as streptolysin O and streptolysin S that can act as antigens. The affected individuals produce specific antibodies against streptolysin O, namely Anti streptolysin O. Determination of these antibodies is very useful for the diagnosis of streptococcal infections and their relative effects such as rheumatic fever and acute glomerulonephritis. An elevated ASO titre of more than 200 IU/ml may indicate an acute streptococcal infection.

## METHOD PRINCIPLE

ASO slide test for detection of antibodies to streptolysin O is based on the principle of agglutination. The test specimen (serum) is mixed with ASO latex reagent and allowed to react. If antibodies to streptolysin O are present then a visible agglutination is observed. If antibodies to streptolysin O are not present or are in concentration less than 200 IU/ml then no agglutination will be observed.

## REAGENTS

Reagent Name	SASO000025T	SASO000050T	SASO000100T
R1 ASO latex	1 vial	1 vial	1 vial
R2 Positive control	1 vial	1 vial	1 vial
R3 Negative control	1 vial	1 vial	1 vial

## WORKING REAGENT PREPARATION AND STABILITY

1. Store the reagent at 2-8°C. DO NOT FREEZE.
2. The shelf life of the reagent is as per the expiry date mentioned on the reagent vial labels.

## SPECIMEN

Only serum should be used for testing. Should a delay in testing occur, store the samples at 2-8°C. Samples can be stored for up to a week. Do not use hemolysed serum.

## MATERIAL PROVIDED WITH THE KIT

**Reagents:** ASO latex reagent, positive control, Negative control

**Accessories:** Glass Slide, Plastic Droppers, Mixing sticks.

## ADDITIONAL MATERIAL REQUIRED

Stop watch, high intensity to room temperature before testing.

## NOTES:

1. In vitro diagnostic reagent for laboratory and professional use only. Not for medicinal use.
2. All the reagents derived from human source have been tested for HBsAg and anti HIV antibodies and are found to be non-reactive. However handle the material as if infectious.
3. Reagent contains 0.1% sodium Azide as preservative. Avoid contact with skin and mucosa. On disposal flush with large quantities of water.

4. The reagent can be damaged due to microbial contamination or on exposure to extreme temperature. It is recommended that the performance of the reagent be verified with the positive and negative controls provided with the kit.
5. The Latex reagent should be shaken well prior to use, to ensure homogeneous suspension of latex.
6. Only a clean and dry glass slide must be used. Clean the slide with distilled water and wipe dry.
7. Accessories provided with the kit only must be used for optimum results.

## PROCEDURE

### QUALITATIVE METHOD

1. Pipette one drop of test sample on to the glass slide using a disposable pipette provided with the kit.
2. Add one drop of ASO latex reagent to the drop of test sample on the slide.
3. Using a mixing stick, mix the serum and the ASO latex reagent uniformly over the entire circle. Do not let the dropper tip touch the liquid on the slide.
4. Immediately start a stopwatch. Rock the slide gently, back and forth, observing for agglutination macroscopically at two minutes. Proceed similarly with each dilution as test specimen.
5. Do not read the results beyond 2 minutes.

### SEMI QUANTITATIVE METHOD

1. Using isotonic saline prepare serial dilutions of the serum sample positive in the qualitative method 1:2, 1:4, 1:8, 1:16 and so on.
2. Pipette the diluted specimens onto the slide. Start with the 1:2 diluted test specimen.
3. Add a drop of ASO reagent to it. Mix well. Spread the mixture uniformly over the entire circle.
4. Immediately start a stop watch. Rock the slide gently, back and forth observing for agglutination macroscopically at two minutes. Proceed similarly with each dilution as test specimen.

## INTERPRETATION OF TEST RESULTS

### Qualitative Method

Agglutination is a positive test result and indicates the presence of detectable levels of Anti-Streptolysin O in the test specimen. No agglutination is a negative test result and indicates the absence of detectable levels of Anti-streptolysin O in the test specimen.

### Semi Quantitative method

Agglutination in the highest serum dilution corresponds to the amount of ASO in IU/ml present in the test specimen.

## CALCULATION

$$ASO \text{ (IU/ml)} = S \times D$$

Where S = Sensitivity of the reagent i.e., 200 IU/ml  
D = Highest dilution of serum showing agglutination.

## REMARKS:

1. Lipemic, hemolysed and contaminated serum samples could produce non-specific results.
2. Serum samples having higher protein content may produce non-specific reagent aggregation.
3. Use of plasma rather than serum can lead to false positive results.
4. Do not read results beyond two minutes.
5. It is recommended that all positive test results should be further tested with methods enabling quantitation of ASO titres.
6. It is recommended that results of the tests should be correlated with clinical findings to arrive at the final diagnosis

## LITERATURE

1. Todd E.W., (1934), J. path and Bact., No.39, 299-320.
2. Klein G.C., (1980) Manual of Clin. Immunol., 7<sup>th</sup> Ed., 431.
3. Spaun J., Bentzon m.w., Larsen S.O. et. Al., (1961) Bull. WHO, 24, 271-279.
4. Klein G.C et., Al., (1971) (Appl. Microbiol., 21.000.

**Genuine Biosystem Private Limited**

Plot No.97 & 98, kattabomman street,  
Parvathy Nagar Extension,  
Old Perungalathur, Chennai - 600063, India.  
Ph: +91-44-48681845

# GB's – CRP

(LATEX AGGLUTINATION TEST)

KIT NAME	KIT SIZE	CAT. NO
GB's CRP	25 Tests	SCRPO00025T
GB's CRP	50 Tests	SCRPO00050T
GB's CRP	100 Tests	SCRPO00100T

## INTRODUCTION

C- Reactive Protein (CRP) is a normal alpha globulin, which increases inflammatory process. The name CRP is derived from the fact that this protein has the capacity to precipitate the somatic C-carbohydrate of Pneumococcus. Elevated CRP levels are usually observed in a variety of infections and inflammatory conditions where there is tissue destruction. The CRP level measurement is useful in differential diagnosis of neonatal septicemia and meningitis. CRP levels are always elevated after myocardial infection and surgery. The CRP test can also help in determining post-surgical complications.

## METHOD PRINCIPLE

Uniform latex particles are coated with anti-human CRP. The specimen containing CRP on mixing with latex Reagent agglutinates, showing a positive test result. If CRP is absent, there will be no agglutination, indicating a negative test result.

## KIT CONTENTS

Reagent Name	SCRPO00025T	SCRPO00050T	SCRPO00100T
R1 CRP latex	1 vial	1 vial	1 vial
R2 Positive control	1 vial	1 vial	1 vial
R3 Negative control	1 vial	1 vial	1 vial

## WORKING REAGENT PREPARATION AND STABILITY

1. Store the reagent at 2-8°C. DO NOT FREEZE.
2. The shelf life of the reagent is as per the expiry date mentioned on the reagent vial labels.

## SPECIMEN

Only serum should be used to testing. Should a delay in testing occur, store the samples at 2-8°C. Samples can be stored for upto a week. Do not use hemolysed serum.

## MATERIAL PROVIDED WITH THE KIT

**Reagents:** CRP latex reagent, positive control, Negative control  
**Accessories:** Glass Slide, Plastic Droppers, Mixing sticks.

## ADDITIONAL MATERIAL REQUIRED

Stop watch, high intensity to room temperature before testing.

## NOTES:

1. In vitro diagnostic reagent for laboratory and professional use only. Not for medicinal use.
2. All the reagents derived from human source have been tested for HBsAg and Anti-HIV antibodies and found to be non-reactive
3. Reagent contains 0.1% sodium Azide as preservative. Avoid contact with skin and mucosa. On disposal flush with large quantities of water.
4. The reagent can be damaged due to microbial contamination or on exposure to extreme temperatures. It is recommended that the performance of the reagent be verified with the positive and negative controls provided with the kit.
5. Shake the CRP latex reagent well before use to disperse the latex particles uniformly and improve test readability.

6. Only a clean and dry glass slide must be used. Clean the slide with distilled water and wipe dry.
7. Accessories provided with the kit only must be used for optimum results.

## PROCEDURE

Bring reagent and samples to room temperature before use.

## QUALITATIVE METHOD

1. Pipette one drop of the test specimen (serum) on the glass slide using disposable pipette provided with the kit.
2. Add one drop of CRP latex reagent to the drop of test specimen on the slide. Do not let the dropper tip touch the liquid on the slide.
3. Using a mixing stick, mix the test specimen and CRP latex reagent uniformly over the entire circle.
4. Immediately start a stop watch, Rock the slide gently back and forth, observing for agglutination macroscopically at two minutes.

## SEMI QUANTITATIVE METHOD

1. Using isotonic saline prepare serial dilutions of the test specimen positive in the qualitative method 1:2, 1:4, 1:8, 1:16, and so on.
2. Pipette each dilution of the test specimen onto separate reaction circles.
3. Add one drop of CRP latex reagent to the drop of test specimen on the slide. Do not let the dropper tip touch the liquid on the slide.
4. Using a mixing stick, mix the test specimen and the latex reagent uniformly over the entire circle.
5. Immediately start a stop watch. Rock the slide gently, back and forth, observing for agglutination macroscopically at two minutes.

## INTERPRETATION OF RESULTS

### Quantitative Method

Agglutination is a positive test result and indicates the presence of detectable levels of CRP in the test specimen.  
 No agglutination is a negative test result and indicates absence of detectable levels of CRP in the test specimen

### Semi Quantitative method;

Agglutination in the highest serum dilution for responds to the amount of CRP in mg/dl present in the specimen.

## CALCULATION

$$\text{CRP (mg/L)} = S \times D$$

Where S = Sensitivity of the reagent i.e., 0.6 mg/L  
 D = Highest dilution of serum showing agglutination.

#### REMARKS:

1. Marked lipemic, hemolysed and contaminated serum samples could produce non-specific results.
2. Use of plasma rather than serum can lead to false positive results.
3. CRP is found to be present after the first trimester of pregnancy and persists until delivery.
4. CRP levels increase in women who are on oral contraceptives.
5. CRP response is not affected by the commonly used anti-inflammatory or immuno suppressive drugs, including steroids, unless the disease activity is affected and it covers an exceptionally broad incremental range up to 3000 times.
6. Do not read results beyond indicated testing time limits.
7. Since CRP production is a non-specific response to tissue injury, it is recommended the result of the test should be correlated with clinical findings to arrive at the final diagnosis.
8. In case where an increase in CRP levels is suspected, but the screening test shows a negative result, semiquantitation should be done to rule out prozone effect.

#### LITERATURE

1. Kidmark, C.O. (1972) Scand J. Clin. Invest. 29, 407.
2. Dey, R.A., Pope, R.M., Perselin, R.H. (1980). J. Rheumatol, 7, 279.

**Genuine Biosystem Private Limited**

Plot No.97 & 98, Kattabomman Street,  
Parvathy Nagar Extension,  
Old Perungalathur, Chennai - 600063, India.  
Ph: +91 44 48681845

## BLOOD GROUPING REAGENT

Ref: 406 BSA 30% 10ml

Ref: 436 BSA 22% 10ml

# Bovine Serum Albumin 22 % and 30% Solution

### Reagent:

Bovine serum albumin (BSA) solutions (22% and 30%) are suitable for the use as a suspension media for antigen determination of human erythrocytes as well as a supplement for compatibility testing, antibody screening and identification. This reagent is finished and intended for laboratory use only.

### Principle of procedure

The procedures used with these reagents are based on the principle of agglutination.

Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed toward the antigen. The BSA solution will act as amplifier medium for existing antibodies.

### Warning

This reagent is prepared from biological material. As biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through exitants of disease. The reagent contains sodium azide, that may be toxic and may react with lead or copper to form highly explosive salts. Because of these reasons reagent should be handled with proper care.

### Storage

Store at 2 to 8°C. May be at room temperature (15 to 30°C) while in use. In principle, store and use the reagents to declared expiry date only.

### Remarks

1. Inappropriate storage impairs efficacy of the reagent.
2. Centrifugation highly different from appointed relative centrifugal force may lead to false results
3. Blood samples to be tested should be within 48 h after collection.
4. Reagents have to be used without additional additives.
5. The procedures identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on automated systems.
6. For usage of this test sera all effective national laws, directives and guidelines have to be observed.

### Reagent preparation

There is no preparation of the reagents required. Use reagents directly from the vials.

### Procedure

Not provided material, additionally needed

#### at Tube Centrifugation Method

1. Test tubes, 10 x 75 mm or 12 x 75 mm
2. Pipettes designed to deliver approximately 100 µL
3. Centrifuge
4. Isotonic saline (with 0,85 - 0,9% sodium chloride)

### Test procedure

#### Tube-centrifugation-test:

(e.g. Antibody identification, screening and compatibility testing)

1. Add one drop of patient's serum in a tube.
2. Add one drop of 2% to 5% suspension of red blood cells in isotonic saline (cells washed one time or up to three times with isotonic saline).
3. Add 2 drops bovine serum albumin solution (22% or 30%) and mix by slightly shaking.
4. Centrifuge for 1 minute at 1.000 rpm (approx. 180 - 270g).
5. Carefully resuspend the mixture and check for agglutination within 3 minutes.
6. Incubate the tube for 15-60 minutes at 37°C.
7. Centrifuge for 1 minute at 1.000 rpm (approx. 180 - 270g).
8. Wash the tube at least 3x with cold isotonic sodium chloride solution.
9. After decanting the tube after the last wash, add 2 drops of anti-human-globulin and mix by slightly shaking.
10. Centrifuge for 1 minute at 1.000 rpm (approx. 180 - 270g).
11. Carefully resuspend the mixture and check for agglutination within 3 minutes.
12. All negative or weak positive tests must be controlled with Coombs control cells.

### Interpretation of results

"Slightly rotating / shaking " at Slide Method and at Tube Centrifugation Method

Positive results (+): visible agglutination of erythrocytes is a positive result and indicates the presence of the corresponding antigen.

Negative results (-): No visible agglutination of erythrocytes is a negative result and indicates the absence of the corresponding antigen.

### Limitations of the test procedure:

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. Bacterial contamination of a component or chemical additive may produce incorrect results.
3. Enzyme treated erythrocytes may react unspecific.
4. If plasma is used instead of serum, complement binding antibodies may not be detected because of missing calcium.

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Monocl. Antisera for Slide and Tube Method

# Anti-A, Anti-B, Anti-AB

## Reagents for ABO Blood Grouping

Murine monoclonal IgM Blood Grouping Reagents.

Store between 2° and 8°C. Inappropriate storage impairs efficacy of the reagent. Preservative < 0,1 % NaN3.

The reagents have been optimised for use as supplied by the recommended techniques without further dilution or additions.

Diagnostic reagents for professional in-vitro diagnostic use only.

## Reagent

The test procedures recommended for this reagent are based upon the agglutination (clumping) of red blood cells which carry a specific antigen in the presence of a corresponding specific antibody. The Cypress-ABO blood grouping reagents contains mouse monoclonal IgM antibodies. When used by the recommended techniques these reagents will cause direct agglutination of red cells carrying the specific antigen. The reagents contain the colours: Anti-A (blue), Anti-B (yellow), Anti-AB (pale yellow).

## Precautions

- The reagent cannot be assumed to be free from infectious agents. Care must be taken in the use and disposal of each container and its contents.
- With each testing, the result of red cells grouping should be confirmed by reverse grouping the individual's serum with known A, and B red cells. Both results have to correspond.
- Centrifugation highly different from appointed relative centrifugal force may lead to false results.
- The procedures identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on automated systems.
- Strength of positive reactions also depends on age of used blood.

## Sample collection

Blood samples should be used as soon as possible. If testing the blood samples is delayed storage should be at 2°-8°C. Samples collected into heparin or oxalate should be tested within two days. Blood drawn into sodium citrate or EDTA should be tested within 14 days. Blood obtained by finger puncture may be tested directly by the slide method but, to avoid clotting, blood collected in this manner should be mixed quickly with the reagent.

## Not provided material, additionally needed:

at Slide Method: Glass slide, Pasteur pipette, Mixing stick

at Tube Centrifugation Method: Test tubes (10x75mm or 12x75 mm), Pipettes designed to deliver approximately 100 µL, Centrifuge, Isotonic saline (with 0,85 - 0,9% sodium chloride)

## Test procedure

### A. Slide Technique

1. Use erythrocyte sediment or whole blood only.
2. On a clean prelabelled glass slide add 1 drop ( $\pm 50 \mu\text{L}$ ) of Cypress ABO reagent and, using a Pasteur pipette, 1 drop ( $\pm 50 \mu\text{L}$ ) of erythrocyte sediment or whole blood.
3. Mix the reagent and the cells with a stick over an area about 2 cm in diameter.
4. By slightly rotating the slide check for agglutination within 1 minute (reaction starts within seconds). Unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.

### B. Tube technique – with centrifugation

1. Use a 2% to 5% suspension of red blood cells in isotonic saline (cells washed one time or up to three times with isotonic saline) only.
2. To 100 µl (alternative 1 drop  $\pm 50 \mu\text{L}$ ) of reagent in a labelled test tube, add an equal volume of suspension of test red blood cells.
3. Mix well by slightly shaking and incubate at room temperature (15 - 30°C) for a 1 – 15 minutes
4. Centrifuge at 1000 rpm ( $\pm 180 - 270 \text{ g}$ ) for 1 minute
5. Gently agitate the tube to dislodge the red cells and examine macroscopically for agglutination within 3 minutes.

## Interpretation of the results

After slightly rotating / shaking at Slide Method and at Tube Centrifugation Method:

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## General Observations

1. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
2. Enzyme treated erythrocytes may react unspecifically.
3. The test serum anti-B doesn't react with "acquired B".
4. Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may give weaker reactivity compared to control cells. Anti-A will not detect all examples of A<sub>x</sub> cells.
5. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test (DAT)) may give weak reactions. In extreme cases, false-negative results may occur.

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## BLOOD GROUPING REAGENTS

vial dropper 10 ml  
Monoclonal IgM + IgG Blended  
Ref 409 Anti D blend  
For Slide and Tube Test

# Anti-D blend

### Intended use

Monoclonal Anti-D blend – testserum is produced from cell culture supernatants of hetero-hybridoma cell lines. One cell line is secreting an antibody of IgG-type and the other antibody of IgM-type that both react specific with the corresponding antigen. The testserum is used to determine whether human red blood cells possess or lack the corresponding blood group antigen D. This testserum reacts negatively with erythrocytes of category D<sup>VI</sup> with the IgM-antibody and positive with the IgG-antibody in indirect Coombs test. The testserum is intended to be used by qualified and technical personnel only. For in vitro use only.

### Principle of procedure

The procedures used with this reagent are based on the principle of agglutination. Normal human erythrocytes, possessing the corresponding antigen, will agglutinate in the presence of the specific antibody directed towards the antigen.

### Reagents

The listed reagent contains antibodies of the following clones:

Anti-D blend (monoclonal, human IgG+IgG, clones: MS-26; TH-28)

The reagent contains < 0,1% (w/v) sodium azide as preservative. Additionally the reagent is prepared of active antibody, sodium chloride, macromolecules and bovine albumin.

### Warning

This reagent was prepared from supernatants of cell cultures. As a biological product it should be looked upon as potentially infectious because of never complete exclusion of danger through exitants of disease. The reagent contains sodium azide, that may be toxic and may react with lead or copper to form highly explosive salts. Because of these reasons this reagent should be handled with proper care.

### Storage

Store at 2 to 8°C. May be at room temperature (15 to 30°C) while in use. In principle, store and use the reagents to indicated expiry date only.

### Remark

1. Strength of positive reactions also depends on age of used blood
2. With each testing positive and negative controls should be performed.
3. Inappropriate storage impairs efficacy of the reagent.
4. Centrifugation highly different from appointed relative centrifugal force may lead to false results.
5. Blood samples to be tested should be used as soon as possible. If a delay in testing occurs, samples should be stored at 2 to 8°C. Blood drawn into heparin or oxalate should be tested within two days. Blood drawn into sodium citrate or EDTA should be tested within 14 days. Blood obtained by finger puncture may be tested directly by the slide method but, to avoid clotting, blood collected in this manner should be mixed quickly with the reagent.
6. The procedures identified below are for manual testing only. When using automated or semi-automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow approved validation procedures to demonstrate compatibility of this product on automated systems.
7. For usage of this testsera all effective national laws, directives and guidelines have to be observed.

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### Procedure

Not provided material, additionally needed:

#### at Slide Method

1. Glass slide
2. Pasteur pipette
3. Mixing stick

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## Test procedure

### Slide Method

1. Use erythrocyte sediment or whole blood only.
2. Place one drop (approximately 50  $\mu\text{L}$ ) of appropriate reagent on a glass slide.
3. Using a Pasteur pipette add one drop erythrocyte sediment or whole blood (approximately 50  $\mu\text{L}$ ) to the glass slide.
4. Mix the erythrocytes with reagent well with a stick and spread to a circle (diameter 2 cm)
5. By slightly rotating the slide, check for agglutination within 1 minute (reaction starts within seconds). Unspecific reactions might appear due to drying of the reaction-formation or if the slide is heated.

### Tube Centrifugation Method

1. Use a 2% to 5% suspension of red blood cells in isotonic saline (cells washed one time or up to three times with isotonic saline) only.
2. Add 100  $\mu\text{L}$  (alternative: one drop = approximately 50  $\mu\text{L}$ ) of appropriate reagent to each tube.
3. Add 100  $\mu\text{L}$  (alternative: one drop = approximately 50  $\mu\text{L}$ ) of appropriate cell suspension to each tube.
4. Mix well by slightly shaking.
5. Incubate tube at room temperature (15 to 30  $^{\circ}\text{C}$ ) for 1 - 15 min.
6. Centrifugation of tube for 1 minute at 2.000 rpm (approximately 800 - 1.000 g).
7. Gently resuspend the red cells and check macroscopically for agglutination within 3 minutes. Document the result.

### Indirect Coombs-test

8. Repeat step 1. - 4. with fresh material or continue after reading of result directly.
9. Incubate tube at 37  $^{\circ}\text{C}$  for 30 min.
10. Wash red cells 3 times with (cold) isotonic saline.
11. Add 100  $\mu\text{L}$  Anti-Human-Globulin-serum (Coombs-serum) to each tube
12. Centrifugation of tube for 1 min at 1.000 rpm (approximately 180 - 270 g).
13. Gently resuspend the red cells and check macroscopically for agglutination within 3 minutes. Document the result.

## Interpretation of the results

" Slightly rotating / shaking " at Slide Method and at Tube Centrifugation Method

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## Limitation of the procedure

1. Inaccuracy at compliance with instructions written under section "Procedures" and "Interpretation of results" may lead to incorrect results.
2. No valid conclusion concerning the test result can be reached, if controls with uncertain or false results occur.
3. Enzyme treated erythrocytes may react unspecifically.
4. Due to variability of antigen expression, reactivity of these reagents against certain phenotypes may give weaker reactivity compared to control cells.
5. Red blood cells coated with alloantibodies or autoantibodies of the same or similar specificity as the reagent (i.e., cells that are positive in the direct antiglobulin-test (DAT)) may give weak reactions. In extreme cases, false-negative results may occur.
6. With slide method most of weak D-antigens ( $\text{D}^{\text{weak}}$ ) and categories are not recognized.
7. This testserum reacts negatively with erythrocytes of category  $\text{D}^{\text{VI}}$  in an agglutinating test and positively in an indirect Coombs-test.
8. The described procedures for testing are for usage with Coombs-serum produced by Cypress Diagnostics. On principle Coombs-sera of other manufacturers can be used, but laboratories have to follow the procedures of this manufacturer. Also laboratories must follow approved validation procedures to demonstrate compatibility of this product.

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1. Use a 2% to 5% suspension of red blood cells in isotonic saline (cells washed one time or up to three times with isotonic saline) only.
2. To 100 µL (alternative 1 drop  $\pm 50 \mu\text{L}$ ) of reagent in a labelled test tube, add an equal volume of suspension of test red blood cells.
3. Mix well by slightly shaking and incubate at room temperature (15 - 30°C) for a 1 - 15 minutes
4. Centrifuge at 1000 rpm ( $\pm 180 - 270 \text{ g}$ ) for 1 minute
4. Gently agitate the tube to dislodge the red cells and examine macroscopically for agglutination within 3 minutes.

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